



Presentation

The cleaning solution PerfectClean is intended for the maintenance of hearing aids and earmolds.

It removes earwax responsible for the decrease in sound quality, increase the life of hearing systems and preserves their effectiveness.

The refill is presented in the form of a container filled with 1.01 fl oz (30mL) of solution, sealed by a tongue and a lid. One refill allows you to complete between 40 and 50 cycles (20 to 25 cycles when the tank is empty during the first use. Quantities given for information only. This data may vary depending on hearing aids).

The patented solution removes completely ear wax. It was developed to be used exclusively with the PerfectClean.

Properties

The specific formulation of the solution complies with the requirements of the bactericidal standard EN 1040 and the yeasticidal standard EN 1275.

The refill meets the requirements of Class I Medical Devices, according to the Regulation 2017/745/EU related to Medical Devices while respecting environmental and health safety requirements.

Using advices

The blue flashing of the PerfectClean indicated that a new refill must be used.

It is recommended to clean your equipment daily.

Once the cycle has ended, rub the hearing aids with a clean absorbent cloth.

Precautionary statements

Read the user manual before use.

Keep out of reach of children.

Product composition

- Water
- Alcohol
- Caprylyl Glycol
- Ethylhexylglycerin
- Benzalkonium Chloride

